

UA-600/600T

Urine Analyzer

- 600 strips per hour
- Cold light source with high luminosity, long shelf life
- Innovative light sensor, high sensitivity
- Automatic detection of strip placing and testing
- Large LCD touch screen or keypad configuration
- Automatic used strip collection
- Standard RS-232 interface for connection with external computer
- Built-in thermal printer and interface for external printer
- Storage of 2000 test results
- Multifunctional software



mindray
healthcare within reach

UA-600/600T

Urine Analyzer

Technical Specifications

| | |
|---------------------|--|
| Throughput | 600 strips per hour |
| Strips | 11 items |
| Test items | Urobilinogen, Bilirubin, Ketone, Blood, Protein, Nitrite, Leukocytes, Glucose, Specific Gravity, Ascorbic Acid, pH |
| Wavelength | Red (R), Green (G) and Blue (B) |
| Principle | Reflectance photometer |
| Display | Large LCD display, Resolution: 240×64 |
| Input | Touch screen (5.7", UA-600T); Keypad (5.3", UA-600) |
| Memory | 2000 results |
| Printer | Built-in thermal printer, paper width: 57mm, interface for external printer |
| Interface | Standard RS-232C |
| Working environment | Temperature 10°C~30°C Humidity ≤85% |
| Power supply | AC 100~240V~, 50/60Hz |
| Power | 45VA |
| Dimension | 390mm(H)×330mm(W)×210mm(D) |
| Weight | 5.5kg |



Accurate Results
Immediate Diagnosis
Convenient

U-11 Urine Strips

Urobilinogen, Bilirubin, Ketone, Blood,
Protein, Nitrite, Leukocytes, Glucose,
Specific Gravity, Ascorbic Acid, pH

High sensitivity and specificity

Result in 6 seconds

Easy operation

UA-600/600T + EH-2050 Plus, help you easily realize the integrated urinalysis



DISTRIBUTOR:



MINDRAY is a trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Specifications subject to changes without prior notice.
© 2011-2012 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.
P/N: ENG-UA600/600T+EH2050Plus-210285/2-20120319

mindray

Mindray is listed on the NYSE under the symbol "MR"
Mindray Building, Keji 12th Road South, High-tech Industrial Park,
Nanshan, Shenzhen 518057, P.R. China
Tel: +86 755 8188 8998 Fax: +86 755 26582680
E-mail: intl-market@mindray.com Website: www.mindray.com