



Rotary Microtome HM 325

Trusted Performance, Added Efficiencies

- High stability and precision for premium-quality paraffin sections
- Unique section thickness setting with removable knobs allowing operation with both the right and left hands
- Resettable electronic section counter
- Two quick-trim stages (10 µm and 30 µm) independent from the preset fine section thickness
- X/Y fine orientation with reproducible zero positioning
- Manual coarse feed wheel
- Specimen retraction (can be deactivated)
- Long stroke length (64 mm) allows high quality sections of Macro / SuperMega cassettes
- Large, removable section waste tray covers the entire working area

Physical Specifications

Dimensions (D x W x H):
19 x 16.5 x 11 in.

(49 x 42 x 28 cm)

Weight: 51 lbs / 23 kg





Horizontal feed range: 28 mm
Section thickness range: 0.5-60 µm
Mechanical trim stages: 10/30 µm
Retraction: 60 µm

Includes

Standard tools, brushed aluminum cover plate, dust cover, and operator guide

Description	Order Number
HM 325 unit only	902100
With Universal Cassette Clamp and E blade holder	902100A
With Universal Cassette Clamp and ER blade carrier	902100ER

You'll need:

-  blades, blade holders, knives
-  Cool Cut and STS
-  microscope slides
-  freezing sprays